

LOW-ASH LIGNITE ANALYSIS USING THE TSI® LIBS DESKTOP ANALYZER

APPLICATION NOTE LIBS-005

Material: Low-ash Lignite

Gippsland Basin (Australia) lignite deposits have characteristically high moisture content (50 to 70%) and a low ash yield (1 to 3%). The lignite is used extensively for domestic electricity production for which 50 Mt p.a is mined.

Test Method

Five crushed and blended lignite samples were prepared and chemically analyzed using acid extraction (ICP-AES analysis). These samples were used to calibrate the TSI® LIBS Desktop Analyzer and then one of the samples was subsequently presented as an unknown five times. Samples for analysis are prepared by placing approximately 5 g of the lignite samples in the sample holder and pressing to 0.3 tone pressure using a LAT hydraulic press.

The analysis time was 20 seconds (all elements) for each sample analyzed.

Detectable Elements

Detectable elements include the principal coal components (C, H, O and N) in addition to the inorganic components Al, Ca, Fe, K, Li, Na, Mg, Si, and Ti.



Detection Limits

Detection limits were determined from three times the standard deviation in multiple measurements of materials of samples with low analyte concentrations. The estimated detection limits for the principal impurities present in low-ash lignite are:

Element	Detection Limit*
Na	0.002 %
Ca	0.002 %
Mg	0.002 %
Fe	0.003 %
K	0.002 %
Ti	0.003 %
Al	0.005 %
Si	0.005 %

**Detection limit expressed in % as-received*

Multiple Analysis Test


All of the following results are expressed as % A.R. (as-received). The moisture content of the test sample was 60% in this case. Some sample heterogeneity was seen in the results with the variances similar to that observed using comparative analysis techniques for this heterogeneous material.

Sample	Al % 309 nm	Ca % 317 nm	Fe % 261 nm	K % 766 nm	Mg % 280 nm	Si % 288 nm	Na % 589 nm	Ti % 325 nm	* S %
1	0.056	0.068	0.060	0.003	0.12	0.11	0.084	0.001	0.39
2	0.048	0.064	0.060	0.003	0.07	0.08	0.072	0.001	0.38
3	0.044	0.060	0.056	0.003	0.09	0.09	0.076	0.0004	0.38
4	0.044	0.052	0.040	0.002	0.08	0.10	0.080	0.001	0.38
5	0.044	0.064	0.064	0.003	0.12	0.08	0.080	0.0004	0.40
Mean (SD)	0.048 (0.004)	0.060 (0.004)	0.056 (0.01)	0.003 (0.0002)	0.10 (0.02)	0.09 (0.01)	0.080 (0.004)	0.001 (0.0005)	0.38 (0.01)
Standard Analysis	0.048 (0.004)	0.064 (0.004)	0.06 (0.002)	0.003 (0.002)	0.08 (0.005)	0.08 (0.005)	0.076 (0.006)	0.002 (0.0002)	0.39 (0.02)

*Sulfur is determined by chemometric correlation to the other detectable components.

Long Term Performance

The TSI LIBS Desktop Analyzer has been installed in commercial lignite-fired power plants and have been operated routinely over a two-year period. Comparison of the TSI LIBS Desktop Analyzer performance to traditional acid-extraction AAS analysis has shown excellent agreement for lignite impurities such as Na, Al, and Ca.



Daily sampling over a one-month period showing the comparison between the concentration of aluminum measured in coal by both AAS and TSI LIBS Desktop Analyzer analysis.



UNDERSTANDING, ACCELERATED

TSI Incorporated – Visit our website www.tsi.com for more information.

USA	Tel: +1 800 874 2811	India	Tel: +91 80 67877200
UK	Tel: +44 149 4 459200	China	Tel: +86 10 8251 6588
France	Tel: +33 4 91 11 87 64	Singapore	Tel: +65 6595 6388
Germany	Tel: +49 241 523030		