



PRODUCT INFORMATION

## Model 1090—Electrical Ionizer

- Charge conditioning to generate a Boltzmann equilibrium charge distribution on aerosol particles without radioactivity.
- Accurate, stable charge conditioning for size distribution measurement by electrical mobility analysis.



### DESCRIPTION

The Model 1090 Electrical Ionizer is an aerosol charge conditioner for conditioning the electrical charge on aerosol particles to a Boltzmann equilibrium charge equivalent to that produced by a radioactive ionizer, but without the use of radioactivity.

The Model 1090 uses an AC corona discharge to generate high concentrations of positive and negative ions for aerosol charge conditioning without generating undesirable contaminant particles. When operated within the design flow rate of 0.5 to 2.5 L/min, the Model 1090 can be used confidently as a replacement for radioactive ionizers for high accuracy aerosol measurement by electrical mobility.

The Model 1090 can be used as a stand-alone ionizer for charge-conditioning in experimental research. It can also be used as a charge-conditioner for aerosol concentration and size distribution measurement by differential mobility and/or scanning mobility analysis. It is also offered as an accessory for the Model 1000XP Wide-Range Particle Spectrometer (WPS) from MSP to measure aerosol concentration and size distribution from 10 to 10,000 nm.

### PERFORMANCE

For aerosols in Boltzmann charge equilibrium, the charge distribution is described by

$$f_n = \frac{\exp(-n^2 e^2 / d_p k T)}{\sum_{n=-\infty}^{\infty} \exp(-n^2 e^2 / d_p k T)}$$

where  $e$  is the elementary unit of charge,  $d_p$  is the particle diameter,  $k$  is the Boltzmann constant,  $T$  is the absolute temperature,  $n$  is the number of elementary units of charge on the particles and  $f_n$  is the fraction of particles in the aerosol carrying  $n$  elementary units of charge. Figure 1 shows the particle charge distribution according to the Boltzmann's law.




Figure 1. Charge distribution according to Boltzmann law

Figure 2 shows the room air size distribution measured by scanning mobility spectrometry using the WPS and the electrical ionizer and Po<sup>210</sup> neutralizer as a charge conditioner.




Figure 3 shows the size distribution analysis of 100.7 nm PSL spheres by scanning mobility spectrometry using the WPS comparing the result using the electrical ionizer and the Po<sup>210</sup> neutralizer.





Figure 4 shows the size distribution of a high concentration NaCl aerosol measured by the WPS when using an electrical ionizer and the Po<sup>210</sup> neutralizer as a charge conditioner.



## APPLICATIONS

- Charge conditioning for electrical mobility aerosol spectrometers
- Aerosol neutralization for aerosol generators
- Aerosol neutralization for filter testing

## SPECIFICATIONS\*

\*Specifications are subject to change without notice

Aerosol Flow Rate	0.5 to 2.5 L/min
Aerosol Diameter	10 nm to 10 µm
Inlet Port	¼" OD Tube
Outlet Port	¼" Swagelok nut with nylon ferrules
Input Voltage	115 or 230 VAC
Input Frequency	50 or 60 Hz
Input Current	<0.2 A @ 115VAC; <0.1 A @ 230 VAC
Dimensions (WxHxD)	25.5 cm x 15.5 cm x 23.2 cm
Weight	5.4 kg.

## TO ORDER

Specify	Description
1090	Electrical Ionizer

### MSP Corporation

5910 Rice Creek Parkway, Suite 300  
Shoreview, Minnesota 55126, U. S. A.  
Phone: 651.287.8100; Fax: 651.287.8140  
sales@mspcorp.com; www.mspcorp.com

Copyright© 2017 MSP Corporation (MSP-PI-1090, Rev. A – A4). All rights reserved.

The MSP logo is a registered trademark of MSP Corporation.  
TSI and the TSI logo are registered trademarks of TSI Incorporated.