



Knowledge Beyond Measure.

FitPro+™ Fit Test Software for PortaCount Pro and Pro+ Respirator Fit Testers



FitPro+™ Fit Test Software makes fit testing simple and foolproof

Included with every PortaCount Pro and PortaCount Pro+ Respirator Fit Tester, the program automatically leads you through the fit test protocol step by step — you can't make a mistake. Your fit test results are stored in a database for easy report generation, retrieval and regulatory compliance.

Computer Requirements

- Microsoft® Windows® 7, Windows 8, Windows 10 or Windows 11 (32-bit & 64-bit) operating systems
- Microsoft .Net version 4 (installation software will install Microsoft .NET automatically, if necessary)
- Note—An Internet connection is required to install Microsoft .NET
- >1 GHz processor
- >1280x800 screen resolution monitor
- 500 MB free space on hard drive
- 1 free USB port
- CD-ROM reader or internet connection (to install the software)

Features and Benefits

- Operate up to four (4) PortaCounts simultaneously using one computer (PortaCount Models 8030/8038 only)
- Select between new Horizontal or Vertical Fit Test Dashboard
- Automates fit test process using built-in fit test exercise protocols
- Microsoft Access® default database (optional support for server-based Microsoft SQL Server®)
- User Interface in five languages (English, German, Polish, French and Spanish)
- Comprehensive database management functionality (merging, moving, copying and deleting records; exporting databases; and backup and restoring of data)
- Maintains separate databases for different departments/locations
- Easily transfers records between PortaCount Fit Testers using a flash drive
- Wide-range of reporting and printing capabilities, including individual fit test reports, due or past due fit test reports, custom reports, and printing of individualized fit test cards (using card-form paper)
- Real-time fit factor display for respirator training and familiarization



Specifications

FitPro+™ Fit Test Software for PortaCount® Pro and Pro+ Respirator Fit Testers

Fit Test Report and Card Printing Options

FitPro+ is designed to improve your respirator program and help protect workers through easy generation of individualized fit test reports inclusive of personal wallet cards. The wallet cards provide documentation of the completed fit test and easy field verification of proper respirator size and model for authorized respirator users.

The reports can be printed by plugging a USB printer directly into the PortaCount Pro/Pro+ Respirator Fit Tester. Exclusive card-form paper is available direct from TSI that provides for the convenient, durable tear-out wallet card.

FIT TEST REPORT

ID NUMBER: 1234
 LAST NAME: BRUCE
 FIRST NAME: MORTON
 COMPANY:
 LOCATION:

TEST DATE: 1/22/2014 10:52 AM
 PRT DATE: 1/22/2015
 PORTACOUNT SN: 80306502
 HW COMPANION: N

RESPIRATOR: MSA ULTRATWIN FULL FACE (900)
 MANUFACTURER: MSA
 MODEL: ULTRATWIN
 MASK STYLE: FULL FACE
 MASK SIZE: F

PROTOCOL: OSHA 29CFR1910.134
 PASS LEVEL: 500
 APPROVAL: EFFICIENCY:95%
 Filter:

EXERCISE	DURATION(S)	FIT FACTOR	PASS
NORMAL BREATHING	64	110391	Y
DEEP BREATHING	64	34756	Y
HEAD SIDE TO SIDE	64	218463	Y
HEAD UP AND DOWN	64	83305	Y
TALKING	64	70531	Y
GRIMACE	15	Excl.	
BENDING OVER	64	165212	Y
NORMAL BREATHING	64	94501	Y
OVERALL FF		77582	Y

FITTEST OPERATOR: _____ DATE: _____
 NAME: JDM
 DATE: _____
 NAME: MORTON BRUCE

NOTE:

Respirator Fit Test Card

Name: MORTON BRUCE Test Date: 1/22/2014
 ID: 1234 Next Test Date: 1/22/2015

Respirator	Results
Mfg: MSA	Overall FF: 77582
Model: ULTRATWIN	FF Pass Level: 500
Style: FULL FACE	Pass: Y
Size: F	Operator: JDM
	Protocol: OSHA 29CFR1910.134
	Fit Test Method: QNFT using TSI PortaCount
*** Your company contact information here ***	

Specifications are subject to change without notice.

Licensed under US Patent No. 5,117,190 and other patents from Commissariat a L'Energie Atomique, France. U.S. Patent No. 6,125,845. Other patents pending.

Access, Microsoft, SQL Server and Windows are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

PortaCount, TSI and the TSI logo are registered trademarks in the United States and may be protected under other country's trademark registrations.



TSI Incorporated - Visit our website www.tsi.com for more information.

USA	Tel: +1 800 874 2811	India	Tel: +91 80 67877200
UK	Tel: +44 149 4 459200	China	Tel: +86 10 8219 7688
France	Tel: +33 1 41 19 21 99	Singapore	Tel: +65 6595 6388
Germany	Tel: +49 241 523030		