

QUICK START GUIDE

HIGH PRESSURE KIT

5300 SERIES

P/N 6011694, REV. B
APRIL 2019




TSI and TSI logo are registered trademarks of
TSI Incorporated.



5300 Series High Pressure Kit (5300-HPK)

Components

Item	TSI Part #	Description	Qty.
1	5300-IF-HP	Inlet Filter—High Pressure	1
2	5300-375NPT	3/8 Inch NPT Tube Ends	2
3	N/A	1/2 Inch Barbed to 3/8 Inch NPT Fitting	1



Assembly

1. Remove collars and replace the existing tube ends with the two (2) supplied 3/8 inch NPT Tube Ends.



2. Screw the 1/2 to 3/8 inch NPT Fitting onto the inlet side of the Filter.



3. Connect outlet side of the Filter to inlet side of 3/8 inch NPT Tube End.



TSI Incorporated – Visit our website www.tsi.com for more information.

USA **Tel:** +1 800 680 1220
UK **Tel:** +44 149 4 459200
France **Tel:** +33 1 41 19 21 99
Germany **Tel:** +49 241 523030

India **Tel:** +91 80 67877200
China **Tel:** +86 10 8219 7688
Singapore **Tel:** +65 6595 6388

P/N 6011694 Rev B

©2019 TSI Incorporated

Printed in U.S.A.

