

MITA MASK ADAPTER KIT MODEL 8125-07

INSTRUCTION GUIDE



APPLICABLE MASK MODELS:
SCOTT AV-2000 & AV-3000 MASKS



Included Accessories

40/74 mm Mask Filter Plug (Qty 1) (for 40 mm threaded (female) adapters) ➤ Installed on mask during Protection Factor test	
Mask Filter Plug (Qty 2) (for ¼-turn adapters) ➤ Installed on mask during Protection Factor test	
Mask Filter Plug (Qty 1) (for 40mm threaded (male) adapters) ➤ Installed on mask during Protection Factor test	

Components Installed on Masks



Example of Scott AV-2000 and AV-3000

From left to right; Example of the 1/4-turn, 40mm threaded (male) and 40mm threaded (female) Scott adapters

The Model 8125-07 MITA Mask Adapter Kit includes all parts necessary for making connections to a variety of Scott filter/cartridge adapters, and is compatible with all Scott AV-2000 and AV-3000 masks which have either a 1/4-turn, 40 mm threaded (male) or 40 mm threaded (female) adapter.



UNDERSTANDING, ACCELERATED

TSI Incorporated – Visit our website www.tsi.com for more information.

USA **Tel:** +1 800 874 2811
UK **Tel:** +44 149 4 459200
France **Tel:** +33 4 91 11 87 64
Germany **Tel:** +49 241 523030

India **Tel:** +91 80 67877200
China **Tel:** +86 10 8219 7688
Singapore **Tel:** +65 6595 6388