

DUSTTRAK™ MOBILE APP

DOWNLOAD DATA WITHOUT
LEAVING YOUR VEHICLE

TSI's Environmental DustTrak™ and various DustTrak models can now communicate with your mobile device through the DustTrak Mobile App. Combined with the wireless router and antenna, the DustTrak Mobile App allows you to quickly download, review and send logged data.


INCREASE PRODUCTIVITY AND
IMPROVE WORKER SAFETY

Outdoor Applications

- + Drive within 150m of the instrument to view and download data.
- + Quickly move on to the next instrument without ever leaving your vehicle.

Indoor Applications

- + View and download data without having to enter the room, reducing your exposure.
- + Data can be accessed through windows, doors and even walls.




UNDERSTANDING, ACCELERATED

SEE WHAT DUSTTRAK MOBILE CAN DO

- + View and graph live mass concentrations
- + Download and review historical data
- + Email data logs directly from the application
- + Delete instrument internal memory, freeing up space for future log files
- + Wireless capability up to 150m with direct line of sight
- + Compatible with Environmental DustTrak and DustTrak models 8530EP, 8533EP, 8540, 8543, 8540-M, 8542-M and 8543-M

FOR MORE INFORMATION, VISIT
tsi.com/DustTrakMobile.
DUSTTRAK MOBILE IS AVAILABLE
FOR DOWNLOAD FROM THE GOOGLE
PLAY™ STORE AND iTunes®.



Immediately view live mass concentrations either numerically or graphically.

Review historical logs prior to sending email.



UNDERSTANDING, ACCELERATED

TSI Incorporated - Visit our website www.tsi.com for more information.

USA Tel: +1 800 874 2811
UK Tel: +44 149 4 459200
France Tel: +33 1 41 19 21 99
Germany Tel: +49 241 523030

India Tel: +91 80 67877200
China Tel: +86 10 8219 7688
Singapore Tel: +65 6595 6388

DustTrak is a trademark, and TSI and the TSI logo are registered trademarks of TSI Incorporated.

Google Play is a trademark of Google Inc.

iTunes is a registered trademark of Apple Inc.