

# MITA Mask Adapter Kit

Model 8125-05



Instruction Guide



Applicable Mask Model:  
MSA Ultra Elite (w/ 40mm threaded, slide- & push-to-connect,  
CBRN versions); MSA Advantage 3000

# Included Accessories

<p><b>40/74 mm Mask Filter Plug</b> (for MSA Ultra Elite masks)</p> <ul style="list-style-type: none"><li>➤ Installed on mask during <b>Protection Factor</b> test</li></ul>	
<p><b>Compression Adapter</b> (for MSA Ultra Elite masks w/ Slide- &amp; Push-To-Connect)</p> <ul style="list-style-type: none"><li>➤ Installed on mask during <b>Protection Factor</b> test</li></ul>	
<p><b>Mask Filter Plug - Qty 2</b> (for MSA Advantage 3000 masks)</p> <ul style="list-style-type: none"><li>➤ Installed on mask during <b>Protection Factor</b> test</li></ul>	

# Components Installed on Masks



40/74 mm  
Mask Filter Plug

## Example of MSA Ultra Elite w/40mm Threaded



40/74 mm  
Mask Filter Plug

Compression  
Adapter

## Example of MSA Ultra Elite w/ Slide- & Push-To-Connect



Mask Filter Plug

### Example of MSA Advantage 3000

TSI and TSI logo are registered trademarks of TSI Incorporated in the United States and may be protected under other country's trademark registrations.



**Knowledge Beyond Measure.**

TSI Incorporated – Visit our website [www.tsi.com](http://www.tsi.com) for more information.

<b>USA</b>	<b>Tel:</b> +1 800 874 2811	<b>India</b>	<b>Tel:</b> +91 80 67877200
<b>UK</b>	<b>Tel:</b> +44 149 4 459200	<b>China</b>	<b>Tel:</b> +86 10 8251 6588
<b>France</b>	<b>Tel:</b> +33 4 91 11 87 64	<b>Singapore</b>	<b>Tel:</b> +65 6595 6388
<b>Germany</b>	<b>Tel:</b> +49 241 523030		

P/N 6006075 Rev. B

©2023 TSI Incorporated

Printed in U.S.A.

