

Standard Accessory Kit

5300 Series



Quick Start Guide

P/N 6011691, Revision C
May 2023



5300 Series Standard Accessory Kit (53000)







Kit Components

TSI Part Number	Description	Picture
5300-05IN	5300 Series 1/2-inch Tube Ends (2)	
5300-IF-LP	5300 Series Inlet Filter Low Pressure	
5000-USBC-A	5000 Series USB-C to USB-A Power and Communications Cable (6 ft.)	
5000-UPS	5000 Series Universal Power Supply	
5000-HUB	5000 Series USB Hub Cable	

5300 Series Optional Accessories

TSI Part Number	Description	Picture
5300-15MM	15-mm ISO Tube Ends (2)	
5300-075IN	3/4-inch Tube Ends (2)	
5300-375NPT	3/8-inch NPT Tube Ends (2) for High Pressure Applications	
5300-EC	5300 Series End Caps (2)	
5300-C	5300 Series Collars (2)	
5300-HPK	5300 Series High Pressure Kit	

5000 Series Optional Accessories

TSI Part Number	Description	Picture
5000-TTS	5000 Series Tabletop Stand	 A tan-colored metal tabletop stand with a flat base and a vertical support. Three screws are shown below it.
5000-RS232	5000 Series USB-RS232 Cable with Null Modem Converter	 A white cable with a USB-A connector on one end and an RS-232 connector on the other. A separate white null modem converter module is shown next to it.
5000-USB-LOCK	5000 Series USB-C to USB-A Cable with Screw Lock (Includes reversible head screwdriver)	 A black USB-C to USB-A cable with a silver screw lock mechanism. A reversible head screwdriver is shown next to it.
5000-WMK	5000 Series Wall Mounting Kit	 A black metal wall mounting kit consisting of a main plate, a smaller plate, and several screws.
5000-WMDRK	5000 Series Wall Mounting DIN Rail Kit	 A black metal wall mounting kit for DIN rail, including a main plate, two DIN rail connectors, and several screws.
5000-DCC	5000 Series Deluxe Carrying Case	 A blue and grey soft carrying case with a handle and a TSI logo on the side.
5000-PC-ADV	5000 Series PC Software (Advanced version)	See following page for more details.

FLO-Sight™ PC Software

Get the most out of your 5000 Series Gas Flow Multi-Meter with TSI® Incorporated's FLO-Sight™ PC Software.

Couple the real-time data visualization of FLO-Sight™ PC software with the fast response rate of our 5000 Series Flow Meters to gain valuable insights and truly understand your flows.

See your flows as you have never seen them before with FLO-Sight software.



FLO-Sight software has these basic features:

- Real-time display of all measurement parameters
- Min, max, and average readings
- Advanced meter configuration
- Volume measurements with user-set begin/end triggers (5310, 5320, 5330 models)

FLO-Sight™ PC Software *(continued)*

FLO-Sight™ software advanced features include:

- Real-time graphing
- Data logging, download, and file management
- Alerts set for any parameter
- Multiple meter operation

A free trial of FLO-Sight™ software's advanced features is included with the Series 5300 Flow Meter. The advanced version can be purchased through the software itself or at www.tsi.com.

Technical Support

For further information/customer support, please contact:

U.S. & INTERNATIONAL

Sales and Customer Service:

(800) 680-1220 / +1 (651) 490-2860

Fax:

+1 (651) 490-3824

TSI Instruments Ltd. (UK)

Sales and Customer Service:

+44 (0) 1494 459200

Email:

answers@tsi.com

TSI and TSI logo are registered trademarks of TSI Incorporated in the United States and may be protected under other country's trademark registrations.



Knowledge Beyond Measure.

TSI Incorporated – Visit our website www.tsi.com for more information.

USA Tel: +1 800 680 1220
UK Tel: +44 149 4 459200
France Tel: +33 1 41 19 21 99
Germany Tel: +49 241 523030

India Tel: +91 80 67877200
China Tel: +86 10 8219 7688
Singapore Tel: +65 6595 6388

